

SMARCD1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 58kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human SMARCD1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

Rsc6p; BAF60A; CRACD1;

Catalog #: A2273

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6602

Isotype: IgG

Swiss Prot: Q96GM5

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the SWI/SNF family of proteins, whose members display helicase and ATPase activities and which are thought to regulate transcription of certain genes by altering the chromatin structure around those genes. The encoded protein is part of the large ATP-dependent chromatin remodeling complex SNF/SWI and has sequence similarity to the yeast Swp73 protein. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).