

EDF1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 16kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human EDF1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

fj

Synonym:

MBF1; EDF-1;

Catalog #: A2283

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8721

Isotype: IgG

Swiss Prot: O60869

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that may regulate endothelial cell differentiation, lipid metabolism, and hormone-induced cardiomyocyte hypertrophy. The encoded protein has also been found to act as a transcriptional coactivator by interconnecting the general transcription factor TATA element-binding protein (TBP) and gene-specific activators. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).