

KAT7

Reactivity: Human Mouse Rat

Tested applications: WB IHC IP

Recommended Dilution: WB 1:200 - 1:2000 IHC 1:20 - 1:200 IP 1:20 - 1:50

Calculated MW: 71kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human KAT7

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

HBO1; HBOA; MYST2; ZC2HC7;

To place an order, please [Click HERE](#).

Catalog #: A2297

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 11143

Isotype: IgG

Swiss Prot: O95251

Purity: Affinity purification

For research use only.