

HINFP

Reactivity: Human Mouse

Tested applications: WB

Calculated MW: 517

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human HINFP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A2306

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 25988

Isotype: IgG

Swiss Prot: Q9BQA5

Purity: Affinity purification

For research use only.