

PRMT8

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 45kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human PRMT8

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

HRMT1L3; HRMT1L4;

Catalog #: A2326

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 56341

Isotype: IgG

Swiss Prot: Q9NR22

Purity: Affinity purification

For research use only.

Background:

Arginine methylation is a widespread posttranslational modification mediated by arginine methyltransferases, such as PRMT8. Arginine methylation is involved in a number of cellular processes, including DNA repair, RNA transcription, signal transduction, protein compartmentalization, and possibly protein translation (Lee et al., 2005)

To place an order, please [Click HERE](#).