

SMYD3

Reactivity: Human

Tested applications: WB FC

Recommended Dilution: WB 1:200 - 1:2000 FC 1:20 - 1:50

Calculated MW: 42kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human SMYD3

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

KMT3E; ZMYND1; ZNFN3A1; bA74P14.1;

Catalog #: A2329

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 64754

Isotype: IgG

Swiss Prot: Q9H7B4

Purity: Affinity purification

For research use only.

Background:

This gene encodes a histone methyltransferase which functions in RNA polymerase II complexes by an interaction with a specific RNA helicase. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).