

E2F1

Reactivity: Human Rat

Tested applications: WB IF IP

Recommended Dilution: WB 1:200 - 1000 IF 1:50 - 1:100 IP 1:20 - 1:50

Calculated MW: 47kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human E2F1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A2442

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 1869

Isotype: IgG

Swiss Prot: Q01094

Purity: Affinity purification

For research use only.