

ATG5

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50- 1:200

Calculated MW: 32kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human ATG5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ASP; APG5; APG5L; hAPG5; APG5-LIKE

To place an order, please [Click HERE](#).

Catalog #: A2468

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 9474

Isotype: IgG

Swiss Prot: Q9H1Y0

Purity: Affinity purification

For research use only.