

EEF1A2

Reactivity: Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:100

Calculated MW: 50kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human EEF1A2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

EEF1A1; EF-1-alpha-2; EF1A; FLJ41696; HS1; STN; STNL;

Catalog #: A2473

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1917

Isotype: IgG

Swiss Prot: Q05639

Purity: Affinity purification

For research use only.

Background:

EEF1A2, also named as EEF1A1, STN and Statin-S1, belongs to the GTP-binding elongation factor family and EF-Tu/EF-1A subfamily. It promotes the GTP-dependent binding of aminoacyl-tRNA to the A-site of ribosomes during protein biosynthesis. EEF1A2 may be critical in the development of ovarian cancer. The MW of EEF1A2 is 50-53kd. This antibody can recognize both EEF1A1 and EEF1A2.

To place an order, please [Click HERE](#).