

GPT2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:100 IF 1:50 - 1:200 IP 1:20 - 1:50

Calculated MW: 58kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human GPT2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ALT2;

Catalog #: A2483

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84706

Isotype: IgG

Swiss Prot: Q8TD30

Purity: Affinity purification

For research use only.

Background:

This gene encodes one of the two pyridoxal enzymes which catalyze a reversible transamination reaction to yield glutamate and pyruvate so these enzymes participate in amino acid metabolism and gluconeogenesis. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2011]

To place an order, please [Click HERE](#).