

RIPK2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 61kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human RIPK2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bp

Synonym:

CCK; RICK; RIP2; CARD3; GIG30; CARDIAK;

Catalog #: A2498

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8767

Isotype: IgG

Swiss Prot: O43353

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the receptor-interacting protein (RIP) family of serine/threonine protein kinases. The encoded protein contains a C-terminal caspase activation and recruitment domain (CARD), and is a component of signaling complexes in both the innate and adaptive immune pathways. It is a potent activator of NF-kappaB and inducer of apoptosis in response to various stimuli.

To place an order, please [Click HERE](#).