

LYZ

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 17kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human LYZ

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LZM;

Catalog #: A2503

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4069

Isotype: IgG

Swiss Prot: P61626

Purity: Affinity purification

For research use only.

Background:

This gene encodes human lysozyme, whose natural substrate is the bacterial cell wall peptidoglycan (cleaving the beta

To place an order, please [Click HERE](#).