

APOBEC3F

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 45kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human APOBEC3F

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

A3F; KA6; ARP8; BK150C2.4.MRNA;

Catalog #: A2507

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 200316

Isotype: IgG

Swiss Prot: Q8IUX4

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the cytidine deaminase gene family. It is one of seven related genes or pseudogenes found in a cluster, thought to result from gene duplication, on chromosome 22.

Members of the cluster encode proteins that are structurally and functionally related to the C to U RNA-editing cytidine deaminase APOBEC1. It is thought that the proteins may be RNA editing enzymes and have roles in growth or cell cycle control. Alternatively spliced transcript variants encoding different isoforms have been identified.

To place an order, please [Click HERE](#).