

HSPB8

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 22kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human HSPB8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

H11; HMN2; CMT2L; DHMN2; E2IG1; HMN2A; HSP22;

Catalog #: A2514

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 26353

Isotype: IgG

Swiss Prot: Q9UJY1

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the superfamily of small heat-shock proteins containing a conservative alpha-crystallin domain at the C-terminal part of the molecule. The expression of this gene is induced by estrogen in estrogen receptor-positive breast cancer cells, and this protein also functions as a chaperone in association with Bag3, a stimulator of macroautophagy. Thus, this gene appears to be involved in regulation of cell proliferation, apoptosis, and carcinogenesis, and mutations in this gene have been associated with different neuromuscular diseases, including Charcot-Marie-Tooth disease.

To place an order, please [Click HERE](#).