

NEDD9

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 92kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human NEDD9

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CAS2; CASL; HEF1; CAS-L; CASS2;

Catalog #: A2521

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4739

Isotype: IgG

Swiss Prot: Q14511

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the CRK-associated substrates family. Members of this family are adhesion docking molecules that mediate protein-protein interactions for signal transduction pathways. This protein is a focal adhesion protein that acts as a scaffold to regulate signaling complexes important in cell attachment, migration and invasion as well as apoptosis and the cell cycle. This protein has also been reported to have a role in cancer metastasis. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).