

GJA4

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:37kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human GJA4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CX37;

Catalog #:A2529

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2701

Isotype:IgG

Swiss Prot:P35212

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the connexin gene family. The encoded protein is a component of gap junctions, which are composed of arrays of intercellular channels that provide a route for the diffusion of low molecular weight materials from cell to cell. Mutations in this gene have been associated with atherosclerosis and a higher risk of myocardial infarction.

To place an order, please [Click HERE](#).