

ADIPOQ

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 26kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human ADIPOQ

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ACDC; ADPN; APM1; APM-1; GBP28; ACRP30; ADIPQTL1;

Catalog #: A2543

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9370

Isotype: IgG

Swiss Prot: Q15848

Purity: Affinity purification

For research use only.

Background:

This gene is expressed in adipose tissue exclusively. It encodes a protein with similarity to collagens X and VIII and complement factor C1q. The encoded protein circulates in the plasma and is involved with metabolic and hormonal processes. Mutations in this gene are associated with adiponectin deficiency. Multiple alternatively spliced variants, encoding the same protein, have been identified.

To place an order, please [Click HERE](#).