

RETN

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 11kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RETN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ADSF; RSTN; XCP1; FIZZ3; RETN1;

Catalog #: A2551

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 56729

Isotype: IgG

Swiss Prot: Q9HD89

Purity: Affinity purification

For research use only.

Background:

This gene belongs to the family defined by the mouse resistin-like genes. The characteristic feature of this family is the C-terminal stretch of 10 cys residues with identical spacing. The mouse homolog of this protein is secreted by adipocytes, and may be the hormone potentially linking obesity to type II diabetes. Alternatively spliced transcript variants encoding the same protein have been found for this gene.

To place an order, please [Click HERE](#).