

IRF2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 39kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human IRF2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

IRF-2;

Catalog #: A2558

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3660

Isotype: IgG

Swiss Prot: P14316

Purity: Affinity purification

For research use only.

Background:

IRF2 encodes interferon regulatory factor 2, a member of the interferon regulatory transcription factor (IRF) family. IRF2 competitively inhibits the IRF1-mediated transcriptional activation of interferons alpha and beta, and presumably other genes that employ IRF1 for transcription activation. However, IRF2 also functions as a transcriptional activator of histone H4.

To place an order, please [Click HERE](#).