

DKK1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 28kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human DKK1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

t

Synonym:

SK; DKK-1;

Catalog #: A2562

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 22943

Isotype: IgG

Swiss Prot: O94907

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that is a member of the dickkopf family. It is a secreted protein with two cysteine rich regions and is involved in embryonic development through its inhibition of the WNT signaling pathway. Elevated levels of DKK1 in bone marrow plasma and peripheral blood is associated with the presence of osteolytic bone lesions in patients with multiple myeloma.

To place an order, please [Click HERE](#).