

PDCD4

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 52kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human PDCD4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

PDCD4; Protein 197/15a; H731; Nuclear antigen H731 like;

Catalog #: A2570

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27250

Isotype: IgG

Swiss Prot: Q53EL6

Purity: Affinity purification

For research use only.

Background:

Programmed cell death 4 (Pdc4) is a novel tumor suppressor that inhibits translation, progression and invasion. It was first identified as being differentially upregulated during apoptosis. Pdc4 interferes with the activity of the eukaryotic initiation factor (eIF) 4A by displacing the scaffold protein eIF4G from its binding to the RNA helicase eIF4A.

To place an order, please [Click HERE](#).