

KRT7

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 51kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human KRT7

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

K7; CK7; SCL; K2C7;

Catalog #: A2574

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3855

Isotype: IgG

Swiss Prot: P08729

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the keratin gene family. The type II cytokeratins consist of basic or neutral proteins which are arranged in pairs of heterotypic keratin chains coexpressed during differentiation of simple and stratified epithelial tissues. This type II cytokeratin is specifically expressed in the simple epithelia lining the cavities of the internal organs and in the gland ducts and blood vessels. The genes encoding the type II cytokeratins are clustered in a region of chromosome 12q12-q13. Alternative splicing may result in several transcript variants; however, not all variants have been fully described.

To place an order, please [Click HERE](#).