

## FES

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 93kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human FES

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

FPS;

**Catalog #:** A2578

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 2242

**Isotype:** IgG

**Swiss Prot:** P07332

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes the human cellular counterpart of a feline sarcoma retrovirus protein with transforming capabilities. The gene product has tyrosine-specific protein kinase activity and that activity is required for maintenance of cellular transformation. Its chromosomal location has linked it to a specific translocation event identified in patients with acute promyelocytic leukemia but it is also involved in normal hematopoiesis as well as growth factor and cytokine receptor signaling. Alternative splicing results in multiple variants encoding different isoforms.

*To place an order, please [Click HERE](#).*