

LAMP1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 44kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human LAMP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

c

Synonym:

LAMPA; CD107a; LGP120;

Catalog #: A2582

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3916

Isotype: IgG

Swiss Prot: P11279

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may also play a role in tumor cell metastasis.

To place an order, please [Click HERE](#).