

BTC

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 20kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human BTC

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A2588

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 685

Isotype: IgG

Swiss Prot: P35070

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the EGF family of growth factors. It is synthesized primarily as a transmembrane precursor, which is then processed to mature molecule by proteolytic events. This protein is a ligand for the EGF receptor.

To place an order, please [Click HERE](#).