

ALDH4A1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 62kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human ALDH4A1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

P5CD; ALDH4; P5CDh

Catalog #: A2595

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8659

Isotype: IgG

Swiss Prot: P30038

Purity: Affinity purification

For research use only.

Background:

This protein belongs to the aldehyde dehydrogenase family of proteins. This enzyme is a mitochondrial matrix NAD-dependent dehydrogenase which catalyzes the second step of the proline degradation pathway, converting pyrroline-5-carboxylate to glutamate. Deficiency of this enzyme is associated with type II hyperprolinemia, an autosomal recessive disorder characterized by accumulation of delta-1-pyrroline-5-carboxylate (P5C) and proline. Alternatively spliced transcript variants encoding different isoforms have been identified for this gene.

To place an order, please [Click HERE](#).