

RHOT2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 68kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RHOT2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RASL; ARHT2; MIRO2; MIRO-2; C16orf39;

Catalog #: A2597

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 89941

Isotype: IgG

Swiss Prot: Q8IXI1

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the Rho family of GTPases. The encoded protein is localized to the outer mitochondrial membrane and plays a role in mitochondrial trafficking and fusion-fission dynamics.

To place an order, please [Click HERE](#).