

NRG4

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 12kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human NRG4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HRG4;

Catalog #: A2599

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 145957

Isotype: IgG

Swiss Prot: Q8WWG1

Purity: Affinity purification

For research use only.

Background:

The neuregulins, including NRG4, activate type-1 growth factor receptors (see EGFR; MIM 131550) to initiating cell-to-cell signaling through tyrosine phosphorylation (Harari et al., 1999 [PubMed 10348342]).

To place an order, please [Click HERE](#).