

LAMTOR5

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 10kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant Protein of human LAMTOR5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LAMTOR5;HBXIP;XIP;Ragulator complex protein LAMTOR5;

Catalog #: A2606

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10542

Isotype: IgG

Swiss Prot: O43504

Purity: Affinity purification

For research use only.

Background:

HBXIP, also named as XIP, belongs to the HBXIP family. It is originally identified by its interaction with the C-terminus of hepatitis B virus (HBV) X protein (HBx). HBXIP is a regulator of centrosome dynamics and cytokinesis. HBXIP is able to up-regulate S100A4 through activating STAT4 and inducing DNA methylation of PTEN. HBXIP is a cellular 18kd protein and cytoplasm stain.

To place an order, please [Click HERE](#).