

MSI1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 39kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human MSI1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MSI1; RNA-binding protein Musashi homolog 1;

Catalog #: A2609

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4440

Isotype: IgG

Swiss Prot: O43347

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein containing two conserved tandem RNA recognition motifs. Similar proteins in other species function as RNA-binding proteins and play central roles in posttranscriptional gene regulation. Expression of this gene has been correlated with the grade of the malignancy and proliferative activity in gliomas and melanomas. A pseudogene for this gene is located on chromosome 11q13.

To place an order, please [Click HERE](#).