

BCOR

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 192kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human BCOR

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MAA2; ANOP2; MCOPS2;

Catalog #: A2646

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 54880

Isotype: IgG

Swiss Prot: Q6W2J9

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene was identified as an interacting corepressor of BCL6, a POZ/zinc finger transcription repressor that is required for germinal center formation and may influence apoptosis. This protein selectively interacts with the POZ domain of BCL6, but not with eight other POZ proteins. Specific class I and II histone deacetylases (HDACs) have been shown to interact with this protein, which suggests a possible link between the two classes of HDACs. Several transcript variants encoding different isoforms have been found for this gene. A pseudogene of this gene is found on chromosome Y.

To place an order, please [Click HERE](#).