

IL11RA

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 45kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human IL11RA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CRSDA;

Catalog #: A2648

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3590

Isotype: IgG

Swiss Prot: Q14626

Purity: Affinity purification

For research use only.

Background:

Interleukin 11 is a stromal cell-derived cytokine that belongs to a family of pleiotropic and redundant cytokines that use the gp130 transducing subunit in their high affinity receptors. This gene encodes the IL-11 receptor, which is a member of the hematopoietic cytokine receptor family. This particular receptor is very similar to ciliary neurotrophic factor, since both contain an extracellular region with a 2-domain structure composed of an immunoglobulin-like domain and a cytokine receptor-like domain. Multiple alternatively spliced transcript variants have been found for this gene.

To place an order, please [Click HERE](#).