

APLNR

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:1000 IHC 1:15 - 1:50

Calculated MW: 43kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human APLNR

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

APJ; APJR; HG11; AGTRL1;

Catalog #: A2659

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 187

Isotype: IgG

Swiss Prot: P35414

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the G protein-coupled receptor gene family. The encoded protein is related to the angiotensin receptor, but is actually an apelin receptor that inhibits adenylate cyclase activity and plays a counter-regulatory role against the pressure action of angiotensin II by exerting hypertensive effect. It functions in the cardiovascular and central nervous systems, in glucose metabolism, in embryonic and tumor angiogenesis and as a human immunodeficiency virus (HIV-1) coreceptor. Two transcript variants resulting from alternative splicing have been identified.

To place an order, please [Click HERE](#).