

KRT15

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 49kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human KRT15

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

K15; CK15; K1CO;

Catalog #: A2660

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3866

Isotype: IgG

Swiss Prot: P19012

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the keratin gene family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. Most of the type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains and are clustered in a region on chromosome 17q21.2. this gene.

To place an order, please [Click HERE](#).