

KRT31

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 47kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human KRT31

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HA1; Ha-1; hHa1; KRTHA1;

Catalog #: A2663

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3881

Isotype: IgG

Swiss Prot: Q15323

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the keratin gene family. As a type I hair keratin, it is an acidic protein which heterodimerizes with type II keratins to form hair and nails. The type I hair keratins are clustered in a region of chromosome 17q12-q21 and have the same direction of transcription.

To place an order, please [Click HERE](#).