

CAP2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:1000 IHC 1:50 - 1:200

Calculated MW: 53kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CAP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A2676

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10486

Isotype: IgG

Swiss Prot: P40123

Purity: Affinity purification

For research use only.

Background:

This gene was identified by its similarity to the gene for human adenylyl cyclase-associated protein. The function of the protein encoded by this gene is unknown. However, the protein appears to be able to interact with adenylyl cyclase-associated protein and actin.

To place an order, please [Click HERE](#).