

CTSZ

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 34kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CTSZ

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CTSX;

Catalog #: A2681

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1522

Isotype: IgG

Swiss Prot: Q9UBR2

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a lysosomal cysteine proteinase and member of the peptidase C1 family. It exhibits both carboxy-monopeptidase and carboxy-dipeptidase activities. The encoded protein has also been known as cathepsin X and cathepsin P. This gene is expressed ubiquitously in cancer cell lines and primary tumors and, like other members of this family, may be involved in tumorigenesis.

To place an order, please [Click HERE](#).