

CCL24

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 13kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human CCL24

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Ckb-6; MPIF2; MPIF-2; SCYA24;

Catalog #: A2684

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6369

Isotype: IgG

Swiss Prot: O00175

Purity: Affinity purification

For research use only.

Background:

This gene belongs to the subfamily of small cytokine CC genes. Cytokines are a family of secreted proteins involved in immunoregulatory and inflammatory processes. The CC cytokines are proteins characterized by two adjacent cysteines. The cytokine encoded by this gene displays chemotactic activity on resting T lymphocytes, a minimal activity on neutrophils, and is negative on monocytes and activated T lymphocytes. The protein is also a strong suppressor of colony formation by a multipotential hematopoietic progenitor cell line. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).