

CDK20

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:100

Calculated MW: 31kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human CDK20

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

P42; CCRK; CDCH; PNQALRE;

Catalog #: A2686

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 23552

Isotype: IgG

Swiss Prot: Q8IZL9

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene contains a kinase domain most closely related to the cyclin-dependent protein kinases. The encoded kinase may activate cyclin-dependent kinase 2 and is involved in cell growth. Alternatively spliced transcript variants encoding distinct isoforms have been reported.

To place an order, please [Click HERE](#).