

CLEC1A

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 32kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CLEC1A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CLEC1; CLEC-1;

Catalog #: A2696

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51267

Isotype: IgG

Swiss Prot: Q8NC01

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the C-type lectin/C-type lectin-like domain (CTL/CTLD) superfamily. Members of this family share a common protein fold and have diverse functions, such as cell adhesion, cell-cell signaling, glycoprotein turnover, and roles in inflammation and immune response. The encoded protein may play a role in regulating dendritic cell function. This gene is closely linked to other CTL/CTLD superfamily members on chromosome 12p13 in the natural killer gene complex region. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).