

CLEC4D

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:1000 - 1:4000 IHC 1:50 - 1:200

Calculated MW: 25kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CLEC4D

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MCL; MPCL; CLEC6; CLEC-6; CLECSF8;

Catalog #: A2697

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 338339

Isotype: IgG

Swiss Prot: Q8WXI8

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the C-type lectin/C-type lectin-like domain (CTL/CTLD) superfamily. Members of this family share a common protein fold and have diverse functions, such as cell adhesion, cell-cell signalling, glycoprotein turnover, and roles in inflammation and immune response. This gene is closely linked to other CTL/CTLD superfamily members on chromosome 12p13 in the natural killer gene complex region. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).