

CNGA2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:10 - 1:50

Calculated MW: 76kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CNGA2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CNCA; CNG2; CNCA1; OCNC1; OCNCa; OCNCALPHA;

Catalog #: A2698

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1260

Isotype: IgG

Swiss Prot: Q16280

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene represents the alpha subunit of a cyclic nucleotide-gated olfactory channel. The encoded protein contains a carboxy-terminal leucine zipper that mediates channel formation.

To place an order, please [Click HERE](#).