

CNTFR

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:400 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 41kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human CNTFR

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CNTFR

Catalog #: A2700

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1271

Isotype: IgG

Swiss Prot: P26992

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the type 1 cytokine receptor family. The encoded protein is the ligand-specific component of a tripartite receptor for ciliary neurotrophic factor, which plays a critical role in neuronal cell survival, differentiation and gene expression. Binding of ciliary neurotrophic factor to the encoded protein recruits the transmembrane components of the receptor, gp130 and leukemia inhibitory factor receptor, facilitating signal transduction. Single nucleotide polymorphisms in this gene may be associated with variations in muscle strength, as well as early onset of eating disorders. Alternatively spliced transcript variants have been observed for this gene.

To place an order, please [Click HERE](#).