

EMC8

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 24kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human EMC8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NOC4; COX4NB; C16orf2; C16orf4; FAM158B;

To place an order, please [Click HERE](#).

Catalog #: A2701

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10328

Isotype: IgG

Swiss Prot: O43402

Purity: Affinity purification

For research use only.