

CDX4

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 30kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CDX4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A2704

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1046

Isotype: IgG

Swiss Prot: O14627

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of a small subfamily of homeobox containing transcription factors. The encoded protein may regulate homeobox gene expression during anteroposterior patterning and hematopoiesis.

To place an order, please [Click HERE](#).