

EEF1G

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 50kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human EEF1G

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

EF1G; GIG35;

Catalog #: A2721

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1937

Isotype: IgG

Swiss Prot: P26641

Purity: Affinity purification

For research use only.

Background:

This gene encodes a subunit of the elongation factor-1 complex, which is responsible for the enzymatic delivery of aminoacyl tRNAs to the ribosome. This subunit contains an N-terminal glutathione transferase domain, which may be involved in regulating the assembly of multisubunit complexes containing this elongation factor and aminoacyl-tRNA synthetases.

To place an order, please [Click HERE](#).