

EGR1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 58kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human EGR1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TIS8; AT225; G0S30; NGFI-A; ZNF225; KROX-24; ZIF-268;

Catalog #: A2722

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1958

Isotype: IgG

Swiss Prot: P18146

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the EGR family of C2H2-type zinc-finger proteins. It is a nuclear protein and functions as a transcriptional regulator. The products of target genes it activates are required for differentiation and mitogenesis. Studies suggest this is a cancer suppressor gene.

To place an order, please [Click HERE](#).