

FOXS1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 35kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human FOXS1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FKHL18; FREAC10;

Catalog #: A2730

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2307

Isotype: IgG

Swiss Prot: O43638

Purity: Affinity purification

For research use only.

Background:

The forkhead family of transcription factors belongs to the winged helix class of DNA-binding proteins. The protein encoded by this intronless gene contains a forkhead domain and is found predominantly in aorta and kidney. The function of the encoded protein is unknown.

To place an order, please [Click HERE](#).