

GNA11

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 42kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human GNA11

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GNA-11;

Catalog #: A2731

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2767

Isotype: IgG

Swiss Prot: P29992

Purity: Affinity purification

For research use only.

Background:

Heterotrimeric guanine nucleotide-binding proteins (G proteins) consist of α , β , and γ subunits and mediate the effects of hormones, neurotransmitters, chemokines and sensory stimuli. To date, over 20 known G subunits have been classified into four families, G(s), G(i/o), G(q) and G(12), based on structural and functional similarities (1,2). Phosphorylation of Tyr356 of G(q)/G(11) is essential for activation of the G protein, since phenylalanine substitution for Tyr356 changes the interaction of G with receptors and abolishes ligand-induced IP3 formation (3).

To place an order, please [Click HERE](#).