

LGALS1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:25 - 1:100

Calculated MW: 15kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human LGALS1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GBP; GAL1;

Catalog #: A2733

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3956

Isotype: IgG

Swiss Prot: P09382

Purity: Affinity purification

For research use only.

Background:

The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. This gene product may act as an autocrine negative growth factor that regulates cell proliferation.

To place an order, please [Click HERE](#).